



CARE DESIGN

Malsch Version:
Item no. 91300 100002.1
DE, Stand 05/2021, Rev. 2.1

Bakare Beds Version:
Rev 1 Initial Issue Updated: 26.4.22



IMPULSE CARE BED

INSTRUCTION MANUAL



Dear Customers,

By purchasing a care bed from Malsch care & clinic design®, you have obtained a long-lasting medical product with functions that meet all the requirements of everyday care while maintaining the highest standards of safety.

Thank you very much for the trust you have placed in us.

Our company guarantees carefully selected materials and continuous quality control while employing state-of-the-art production technologies.

Complying with the usage and operating instructions helps to prevent the risk of accidents and preserves the high value of your care bed.

CONTENTS

Product range	6
Specific function	7
Environmental sustainability	7
Note regarding the instructions manual	8
Pictograms / symbols	8
Safety instructions	9
Technical data	14
Rating plate	16
Hand controller symbols	19
IMPULSE care bed with undercarriage Edition 100 / Edition 400 / Edition 400 LR / Edition XL (design <u>with</u> axial tilting sleeping surface)	19
IMPULSE care bed with undercarriage Edition XL (design <u>without</u> axial tilting sleeping surface)	20
IMPULSE care bed with undercarriage Edition 420 / Edition 420 LR	21

Functions	22
IMPULSE care bed with undercarriage Edition 400	23
IMPULSE care bed with undercarriage Edition 400 LR	24
IMPULSE care bed with undercarriage Edition 420	25
IMPULSE care bed with undercarriage Edition 420 LR	26
IMPULSE care bed with undercarriage Edition XL	27
Description of function	28
Back rest	28
Mechanical back rest release / CPR (option)	28
Upper / lower leg section	29
Auto contour	29
Height adjustment	29
Transfer position	29
Low position / fall prevention	30

Trendelenburg position	30	Disposal instructions	48
Anti-Trendelenburg position	31	Accessories (optional)	49
Comfort sitting position	31	Trapeze bar	49
Hand controller locking function	31	IV drip holder	49
Braking and transport	32	Bedside light	49
Adjust DS side rails	34	Hand controller bracket	50
Adjust GS side rails	36	Integrated bed linen holder	50
Adjust VGS side rails	38	Troubleshooting	51
Adjust EVGS side rails	39	Product safety	52
VGS side rail dimensions	40	Cleaning and disinfection	53
Integrated bed extension	42	Safe decommissioning/storage	54
Maintenance	44	Electromagnetic compatibility (EMC)	55
Procedure	44	Warranty and service	58
Battery replacement	46	Customer service	58
Maintenance intervals	47	Declaration of conformity	58
Delivery and assembly	47	Certificates	59

PRODUCT RANGE



IMPULSE care bed with undercarriage
Edition 100 "Style"



IMPULSE care bed with undercarriage
Edition XL



IMPULSE care bed with undercarriage
Edition 400



IMPULSE care bed with undercarriage
Edition 400 LR



IMPULSE care bed with undercarriage
Edition 420



IMPULSE care bed with undercarriage
Edition 420 LR

SPECIFIC FUNCTION

The care beds made by Betten Malsch GmbH are used in the care of residents with physical limitations. The beds are designed exclusively for this purpose. The functions of these care beds assist the care staff in their daily work and offer convenient solutions for positioning the resident and compensating for certain symptoms experienced by residents of retirement and care homes or comparable medical facilities. This complies with working environment 3 and 5 as stipulated by IEC 60601-2-52:2009/AMD1:2015.

Prior written consent from Betten Malsch GmbH is required if the care beds are to be used for other applications.

The product is intended for use as a care aid or health device. As such, it is subject to the regulations of the relevant insurance associations. The care bed is a medical product with reference to applicable industry standards and regulations. Therefore, this product must only be used under medical supervision.

The care beds described in this instruction manual are intended for adult residents with a body weight of at least 40 kg and a height of at least 146 cm. In accordance with the standard IEC 60601-2-52:2009/AMD1:2015, the beds must not be used by residents whose body weight and height are below these limits or who have a BMI under 17, as the risk of injury is significantly increased for this group.



Caution! The use of incompatible side rails and mattresses can lead to injury as body parts may become trapped.

ENVIRONMENTAL SUSTAINABILITY

Betten Malsch GmbH care beds are manufactured in line with the relevant regulations using state-of-the-art processing technologies and are free from hazardous materials. The materials used to finish surfaces are CFC- and solvent-free.

Care beds that are taken out of service due to their age or irreparable damage must be disposed of in line with local disposal regulations.



Caution! Observe the relevant local regulations when disposing of metal, wood and electrical waste.

NOTE REGARDING THE INSTRUCTIONS MANUAL

The following instructions and guidelines in this instruction manual are intended for care staff or other persons and staff tasked with operating and using the care bed.



The instruction manual must be accessible to personnel at all times to avoid operating errors and to guarantee fault-free operation. The care staff must have a good understanding of the care beds and be trained in their operation before using them for the first time. The instruction manual must be used for this training.

The instruction manual has been written for the IMPULSE care bed. The images, graphics and texts they contain may differ from the equipment supplied.



The manufacturer offers technician training for maintenance and servicing work on their care beds. A certificate obtained as part of this training authorises the holder to carry out technical work independently on the beds.

PICTOGRAMS / SYMBOLS

For better orientation, this instruction manual uses the pictograms described below.



Important!

Instructions labelled in this way must be strictly observed in order to avoid injury or damage!



Information!

This pictogram identifies information relating to the current subject.

SAFETY INSTRUCTIONS

It is important that the following safety instructions are observed to prevent risks to residents and carers and to avoid any damage to the bed:

- ⚠ The instruction manual must be read and observed before using the care bed.
- ⚠ It is vital to observe the information given on the rating plate! The information on the rating plate is explained in detail on *p. 16* of this instruction manual.
- ⚠ In the event of any faults or defects that could endanger persons, the bed must not be used.
- ⚠ Electrically adjustable care beds must only be operated by the resident after instruction by trained staff.
- ⚠ Before the bed is used for the first time, the operator must ensure that it is safe to use and in good condition.
- ⚠ The castors must always be set to the braked position to ensure the resident does not fall when getting into or out of the bed.
- ⚠ The bed can be moved into various positions. When doing so, take care to ensure no parts of the body or other objects are located in the adjustment area.
- ⚠ Only care staff are permitted to adjust the side rails. When adjusting the position of the sleeping surface, take care to ensure the resident does not come into contact with the side rails to avoid trapping any part of the body.

- ⚠ The functionality of the side rails must be checked every day. They must not bear any load of over 75 kg vertically or over 50 kg horizontally.
- ⚠ When using CPR (optional, mechanical emergency lowering of the back rest), always additionally relieve the load on the back rest by hand to prevent the back rest from dropping in an uncontrolled manner.
- ⚠ The IMPULSE care bed features a battery-operated emergency mode. This allows the one-off lowering of the sleeping surface in the event of a power cut. The 9-V batteries in the control unit must be checked during the annual safety inspection and replaced at least every two years.
- ⚠ The hand controller functions can be locked or released on the rear side using the key switch. Check that the locking function has taken effect on the hand controller.
p. 19–20 Hand controller symbols
- ⚠ The drive system used must be – operated using a VDE-approved power source – 100–240 V, 50/60 Hz mains socket.
- ⚠ The mains connection cable is also protected by a mechanical strain relief device. Nevertheless, take care to ensure that no sharp edges, mechanical stresses or pinch/shear points are present.
- ⚠ The design of the hand controller means it can be hung on the bed in such a way that the buttons are not activated inadvertently by squashing between two objects. Take care to ensure that the hand controller is fully accessible and not trapped between the side rail and bedside table.

⚠ Observe safety distances to walls, window ledges and other furnishings when using the care bed in an resident's room. The safety distances depend on the design and model of the care bed and are based on the height adjustment and the tilting motions of the bed. The minimum distance is 30 mm.

⚠ Improper use of the bed may cause hazards. Examples of improper use include:

- Unauthorised activation of the electrical functions
- Use of the bed by persons with a body weight of less than 40 kg or a BMI of less than 17, or a height of less than 146 cm *☞ p. 7 Specific function*
- Moving the bed by pulling on the mains cable or side rails
- More than one person adjusting the bed at the same time
- Activation of the functions by the resident without prior instruction
- Pulling the mains cable to disconnect it from the power supply
- Moving the bed on sloping or unsurfaced ground

⚠ In line with the standard DIN EN 60601-2-52:2016, when choosing a mattress, it is important that there is a minimum distance of 22 cm between the top of the sleeping surface and the top of the side rails in its fully extended position. The mattress used must meet the applicable safety standards.

⚠ The constant presence of liquid in the area of the motor must be avoided (e.g. incontinence).

⚠ For safety reasons, the handle on the trapeze bar must be replaced completely every 5 years.

⚠ Servicing and repairs on electrical components must be carried out by specially trained staff and only original replacement parts from the manufacturer must be used.

⚠ The care bed is not suitable for extended operation beyond a working cycle of two minutes. If the mains adapter is overloaded or if it overheats, it will shut off automatically. Further operation is possible only after a 30-minute cool-down phase. Observe the drive manufacturer's notes on the rating plate!

⚠ It is essential to avoid obstructing any part of the bed mechanism as this can lead to damage or complete disabling of the drive system due to overheating.

⚠ Likewise, the safe working load must not be exceeded.

⚠ If an immobile resident remains in the same position for a long period without the use of additional positioning aids, this can lead to pressure sores. The manufacturer of the care bed is not liable for this in any way.

⚠ Electrically operated care beds are active medical products and must be maintained according to Article 7 of the Medical Devices Operator Ordinance (MedProd-BetrV). These maintenance measures must be carried out regularly (at

least once per year). This must involve visual and operational inspections of functional and electrical safety in line with VDE0751. *☞ p. 46 Maintenance*

⚠ Furthermore, electrically operated care beds are electrical appliances and their safety is the responsibility of the employer. The supervisory function of this obligation is the responsibility of the Employers' Liability Insurance Association for Health Service and Welfare Work (BGW) and the Trade Supervisory Board (Gewerbeaufsichtsamt). The regulations of the employers' liability insurance associations apply, particularly those of the German statutory accident insurance body (DGUV), rule 3 of which stipulates regular inspections of movable electrical equipment at a recommended interval of six months, but at least once a year. These inspections may only be carried out by a certified electrician or person with electrical training using specialist measurement and inspection equipment. The inspections according to DGUV rule 3 can be conducted by specialist staff trained by the manufacturer as part of the inspections and maintenance service for medical products.

⚠ Electrically operated care beds are active medical products and must be listed in an inventory for each site in line with Article 13 of the German Medical Devices Operator Ordinance (MedProd-BetrV). It is advisable to also document the correct implementation of the required technical safety checks in this inventory and to specify the date of the next inspection. The required protocols concerning technical safety checks already performed must be appended to the inventory.

⚠ Proper execution and traceable documentation of the technical checks, maintenance and servicing work prescribed by the manufacturer, as well as the technical safety checks, are required in order to preserve the warranty rights of the purchaser. If the operator of a medical product does not meet their obligations, this could lead to the risk of damage and accidents for which the manufacturer is explicitly not liable.

⚠ Maintenance work must be carried out and documented by trained staff.

⚠ The bed must be left in the lowest position if the resident is unattended to reduce the risk of injury caused by falling out of bed.

⚠ If the mains connection cable is damaged, the bed must no longer be used and must immediately be taken out of operation.

⚠ Improper use of the mains connection cable can result in hazards (e.g. electric shock). Examples of this are cable breaks due to kinking, shearing or other mechanical damage.

⚠ When using other ME devices in conjunction with the care bed, precautions must be taken to prevent damage to the device cable or other components of the ME device resulting from crushing between the movable parts of the medically used bed.

⚠ The care bed is not suitable for use in the vicinity of active facilities that use high-frequency surgical devices.

⚠ The care bed is not suitable for use in HF-shielded rooms used for magnetic resonance imaging in which high-intensity EM disturbance variables occur.

⚠ The use of this care bed directly next to or in conjunction with other electrical devices (e.g. stacked) must be avoided as this can lead to faulty operation. If using the bed in the manner described above is absolutely necessary, the devices involved should be subjected to a function test for a longer period of time, in order to rule out malfunction due to interference.

⚠ The use of other accessories, transducers and cables other than those specified or provided by Betten Malsch GmbH can cause increased emissions of electromagnetic interference or reduce the electromagnetic interference resistance of the device and lead to faulty operation.

⚠ Portable HF communication devices (radio devices) – including their accessories, such as antenna cables and external antennas – must not be used closer than 30cm to the parts and cables of care beds as designated by Betten Malsch GmbH. Non-observance can impair the performance of the care bed.

⚠ The emissions of this device are below the thresholds defined by IEC/CISPR 11:2009, Class A and permit the use of it in industrial environments and hospitals. This device may not provide adequate protection against wireless services if used in residential areas (for which class B is normally

required according to CISPR 11). The user may have to implement remedial measures such as relocating or repositioning the device.

⚠ Servicing and maintenance tasks may not be performed while the ME device is in use.

TECHNICAL DATA

Model	Dimensions [cm]	Sleeping surface [cm]	Height adjustment [cm]	Weight	Maximum load	Anti- / Trendelenburg ²	Upper leg adjustment	Back rest adjustment
IMPULSE care bed with undercarriage Edition 100	206 x 100	200 x 90	37 to 84°	approx. 120 kg	225 kg total 190 kg resident 20 kg mattress 15 kg accessories	17° / 14°	33°	71° 12 cm mattress compensation
	206 x 110	200 x 100						
IMPULSE care bed with undercarriage Edition 400 Edition 400 ZB	206 x 100	200 x 90	25 to 82°	approx. 120 kg	225 kg total 190 kg resident 20 kg mattress 15 kg accessories	17° / 14°	33°	71° 12 cm mattress compensation
	174.5 x 90	168.5 x 80						
	186 x 90	180 x 80						
	196 x 90	190 x 80						
	196 x 100	190 x 90						
	206 x 90	200 x 80						
	206 x 120	200 x 110						
	206 x 130	200 x 120						
	226 x 100	220 x 90						
IMPULSE care bed with undercarriage Edition 400 LR	206 x 100	200 x 90	26 to 82°	approx. 120 kg	225 kg total 190 kg resident 20 kg mattress 15 kg accessories	17° / 14°	33°	71° 12 cm mattress compensation
IMPULSE care bed with undercarriage Edition 420	206 x 100	200 x 90	27 to 80°	approx. 100 kg	200 kg safe working load 165 kg resident	–	30°	70° 12 cm mattress compensation
	206 x 90	200 x 80						
IMPULSE care bed with undercarriage Edition 420 LR	206 x 100	200 x 90	28 to 81°	approx. 84 kg	200 kg safe working load 165 kg resident	–	30°	70° 12 cm mattress compensation
	206 x 90	200 x 80						
IMPULSE care bed with undercarriage Edition XL	206 x 110	200 x 100	33 to 83°	approx. 160 kg	300 kg total 250 kg resident 25 kg mattress 25 kg accessories	17° / 14°	33°	71° 12 cm mattress compensation
	226 x 110	220 x 100						

*measured from sleeping surface frame

LOW-VOLTAGE DRIVE SYSTEM

(SMPS external switching power supply)

Model	IMPULSE care bed with undercarriage Edition 100	IMPULSE care bed with undercarriage Edition 400 Edition 400 ZB	IMPULSE care bed with undercarriage Edition 400 LR	IMPULSE care bed with undercarriage Edition 420 Edition 420 LR	IMPULSE care bed with undercarriage Edition XL
Electrical connection	100 - 240 V AC 50 / 60 Hz				
Output voltage	35 V DC 2 A				35 V DC 2,5 A
Over-current off	7.5 - 11.5 A DC				
Over-voltage off	45 V DC				
Standby operation	max. 0,5 W				
Protection	IPX 4				
Protection class	II				


Pressure force

Lifting system force	2 x 3000 N	2 x 3000 N	1 x 8000 N	2 x 6000 N
RS adjustment force	2 x 3000 N	2 x 3000 N	2 x 4000 N	3 x 3000 N
Motor running time	2 min. / off 18 min.			

Data on operation, transport and storage


Operating temperature range	+10° C to +40° C
Transport/storage temperature range	+5° C to +50° C
Relative humidity	30 % to 75 %
Atmospheric pressure range	700 hPa to 1060 hPa
Operating volume	54 dB (A)
Operating altitude	max. 3000 m


RATING PLATE




1  **Betten Malsch GmbH | Rohbergstraße 9 | D-36208 Wildeck-Obersuhl**
2 Tel.: +49 (0)6626 915-100 | Fax: +49 (0)6626 915-116








3 **REF** **Care bed IMPULSE**
SN **0819 1234567**

4 **Input:** 100-240V AC 50/60 Hz 2.1-0.9 A
Output: 35 V DC 2/2.5 A
5 **Operation:** max. T_on: 2 min.
min. T_off: 18 min.
6 **Protection class:** IPX4





7 


9  = 225 kg
 = 190 kg
 = 375 kg

8       




Example illustration of a rating plate of an IMPULSE without UDI








1  **Betten Malsch GmbH | Rohbergstraße 9 | D-36208 Wildeck-Obersuhl**
2 Tel.: +49 (0) 6626 915-100 | info@bettenmalsch.de | bettenmalsch.com

3 **REF** **Care bed IMPULSE**
SN **0121 0400 1234567**


10 **UDI**  **11**
(xx)xxxxxxxxxxxxxx
(xxx)xxxxxxxxxxxxxx

4 **Input:** 100-240 V AC 50 / 60 Hz 2.1 - 0.9 A
5 **Output:** 35 V DC 2 A
6 **Operation:** max. 2 min. ON / 18 min. OFF
7 **Protection class:** IPX4

9  = 225 kg
 = 190 kg
 = 375 kg



8       






Example illustration of a rating plate of an IMPULSE with UDI

The rating plate is located on the underside of the sleeping surface, at the head of the bed on the right. To inspect the rating plate, raise the head section to the upper position.




1. Address of manufacturer
2. Model ID
3. Serial number
4. Electrical voltage, frequency, power consumption
5. Operating time of the motorised adjustment system: please take note of this information to prevent overheating. For example, the bed must not be operated continuously for more than 2 minutes, and after 2 minutes of continuous operation, it must be switched off for at least 18 minutes.
6. Protection of electrical equipment against splash water "Only use in dry areas"; protection rating class II (double insulation, protective insulation)
7. Indicates the next technical check after delivery in line with VDE0751-1

8. Explanation of the symbols used on the rating plate:

	Labelling as medical product
	Application part type B

	Directive 2002/96/EC relating to old electrical and electronic equipment
	Conformity marking in line with the Medical Device Directive
	Protection rating class IPX4 and protection rating class II
	"Only use in dry areas"
	"Observe the instruction manual"

9. Explanation of the symbols used on the rating plate:

	Safe working load
	Maximum permissible weight of the resident
	Maximum total weight of the medical product including resident (Bed weight plus maximum permissible weight of resident)


10. 2D-Barcode (GS1 Data Matrix) DI + PI = UDI

11. DI (Device Identifier)
PI (Product Information Number)

Labels


A separate sticker to the right of the rating plate refers to the labels described below:

Symbol: Label indicating beds for adults used for medical purposes in line with IEC 60601-2-52:2009/AMD:2015



The care bed is approved for adult residents with a body weight of at least 40kg and a height of at least 146cm. In accordance with the standard IEC 60601-2-52:2009/AMD1:2015, the care bed must not be used by residents whose body weight and height are below these limits or who have a BMI under 17, as the risk of injury is significantly increased for this group.

Symbol: Label indicating replaceable mattresses in line with IEC 60601-2-52:2009/AMD:2015 – please observe the information and instruction manuals for the mattresses.




The following table contains information concerning mattress sizes depending on the sleeping surface dimensions:

Mattress size [cm]	Sleeping surface dimensions [cm]	Volumetric weight [kg/m ³]
78×200×12/14	80×200*	35 - 50
88×200×12/14	90×200	35 - 50
98×200×12/14	100×200*	35 - 50
108×200×12/14	110×200*	35 - 50
118×200×12/14	120×200*	35 - 50

* Optional special sizes


HAND CONTROLLER SYMBOLS HC-146

IMPULSE care bed with undercarriage Edition 100/ Edition 400/ Edition 400LR/ Edition XL (design with axial tilting sleeping surface)




Locking function

Positioned on the back of the hand controller for restricting operation by resident.




HAND CONTROLLER SYMBOLS HC-146

IMPULSE care bed with undercarriage Edition XL
(design without axial tilting sleeping surface)




Locking function

Positioned on the back of the hand controller for restricting operation by resident.




HAND CONTROLLER SYMBOLS HC-147

IMPULSE care bed with undercarriage
Edition 420/Edition 420LR




Locking function

Positioned on the back of the hand controller for restricting operation by resident.



HAND CONTROLLER SYMBOLS HB-400

IMPULSE care bed with undercarriage Editionen 400/400LR/500



Locking function


To activate the locking function, place the magnetic chip supplied in the area marked red. At the same time, press the respective blocking button (for secondary symbols, see table). A red lighted LED in the middle of the row of buttons confirms that the function has been locked.

Button	Locking
S2	all functions
S10	Trendelenburg position
S12	low position

The locked functions are released in the same way. By placing the magnetic chip and simultaneously pressing the button with the lock symbol, the function is released again and the respective LED goes out.

FUNCTIONS

IMPULSE care bed with undercarriage Edition 100




1. Head/foot board
2. Adjustable lower leg section
3. Adjustable upper leg section
4. Seat section (fixed)
5. Adjustable back rest
6. Trapeze bar mounts on both sides at head end of sleeping surface
7. Stand with electronic castor positioning

FUNCTIONS

IMPULSE care bed with undercarriage Edition 400

Design with vertically telescopic side rails with one-handed operation, split on both sides (EVGS 7.5.31)




1. Head/foot board
2. Integrated two-stage bed extension (+ 10/+ 20 cm)
3. Adjustable lower leg section
4. Adjustable upper leg section
5. Seat section (fixed)
6. Adjustable back rest
7. Trapeze bar mounts on both sides at head and foot end of sleeping surface
8. Mechanical CPR emergency lowering of the back rest (optional)
9. Telescopic side rails (VGS)
10. Unlocking side rails (one-handed operation)
11. Telescopic side rail movement
12. Lower leg section adjustable notched bracket
13. Foot pedal for arresting each axle (▲ ▼) of covered 50mm castors

FUNCTIONS

IMPULSE care bed with undercarriage Edition 400LR


Design with vertically telescopic side rails with one-handed operation, head-end (EVGSK 7.5.31)



FUNCTIONS

IMPULSE care bed with undercarriage Edition 420

Design with swivelling side rails, split on both sides (GS V3)




1. Head/foot board
2. Integrated two-stage bed extension (+ 10/+ 20 cm)
3. Adjustable lower leg section
4. Adjustable upper leg section
5. Seat section (fixed)
6. Adjustable back rest
7. Trapeze bar mounts on both sides at head and foot end of sleeping surface
8. Swivelling side rails (GS V2¹)
9. Telescopic side rail movement
10. Lower leg section adjustable notched bracket
11. Foot pedal for arresting each axle (▲ ▼) of covered 50 mm castors

¹ Only on versions without integrated bed extension, otherwise V3

FUNCTIONS

IMPULSE care bed with undercarriage Edition 420 LR

Design with swivelling side rails, split on both sides (GS V3)




1. Head/foot board
2. Integrated two-stage bed extension (+ 10/ + 20 cm)
3. Adjustable lower leg section
4. Adjustable upper leg section
5. Seat section (fixed)
6. Adjustable back rest
7. Trapeze bar mounts on both sides at head and foot end of sleeping surface
8. Swivelling side rails (GS V2¹)
9. Telescopic side rail movement
10. 125 mm double castors
11. Lower leg section adjustable notched bracket
12. Foot pedal with central brake (▼), fixed direction (▲) and free wheel (→) of the double castors

¹ Only on variants without integrated bed extension, otherwise V3

FUNCTIONS

IMPULSE care bed with undercarriage Edition XL

Design with swivelling side rails, split on both sides (GS V3)



1. Head/foot board
2. Integrated two-stage bed extension (+ 10/+ 20 cm)
3. Adjustable lower leg section
4. Adjustable upper leg section
5. Seat section (fixed)
6. Adjustable back rest
7. Trapeze bar mounts on both sides at head and foot end of sleeping surface
8. Swivelling side rails (GS V2¹)
9. Telescopic side rail movement
10. CPR emergency lowering (optional, design with GS)
11. Lower leg section adjustable notched bracket
12. 125 mm double castors
13. Foot pedal with central brake (▼), fixed direction (▲) and free wheel (—) of the double castors

¹ Only on variants without integrated bed extension, otherwise V3

DESCRIPTION OF FUNCTION

Depending on the design and type of care bed, the bed can be adjusted into various positions by adjusting the back rest, the upper leg section, the lower leg section and the height of the bed:

Back rest

Use the corresponding buttons on the hand controller to adjust the back rest.



(Back rest operating buttons)

When raised, the back rest can be adjusted by 120mm over its usual length to the head board or 100mm with Edt. 420/420 LR (optional). Undercarriage Edition 420/420 LR optionally available with 100mm mattress compensation of the backrest.

This function (mattress compensation) allows the resident to sit in a comfortable position while avoiding compression of the buttocks area and without compressing or restricting the stomach or upper body.

Mechanical back rest release/CPR (option)

Model for VGS

Handle at head end on a level with the side rails.

When operating the mechanical release, hold the backrest firmly in the current



position and, if possible, relieve some of the load. Pull the release to disengage the back rest and lower it manually into the end position.

Model for GS/ Model without side rail

Handle at the middle of the bed underneath the side panel.

When operating the mechanical release, hold the backrest firmly in the current position and, if possible, relieve some of the load. Pull the release to disengage the back rest and lower it manually into the end position.



Press the S1 button on the hand controller again to reactivate electrical back rest adjustment.



Caution! Before pressing the lever, ensure that there are no objects or parts of the body below the back rest. Manually relieve the weight on the back rest during adjustment to prevent it falling in an uncontrolled manner.

¹ Undercarriage Edition 400/400LR/100/XL

Upper/lower leg section

Use the corresponding buttons on the hand controller to adjust the upper leg section.



Care staff can move the lower leg section into a horizontal position (extended leg elevation) using the adjustable notched bracket.

Auto contour

Use the corresponding buttons on the hand controller to adjust the auto contour.



Using the button function adjusts the back rest and upper leg section equally.

This function must only be used with mobile residents and residents without any physical problems.



Caution! Please consider resident safety when adjusting the auto contour! Ensure that no body parts or any objects are located in the area of the lifting mechanism.

Height adjustment

Use the corresponding buttons on the hand controller to adjust the height.



Caution! Please consider resident safety when adjusting the height. Ensure that no body parts or any objects are located in the area of the lifting mechanism.

Transfer position

(only with IMPULSE care bed with undercarriage Edition 400 and IMPULSE care bed with undercarriage Edition 400 LR)

The transfer position facilitates optimum mounting and dismounting with the seat area of the bed at a sitting height. Hold the operating button down until the position is reached.



Low position/fall prevention

(only with IMPULSE care bed with undercarriage Edition 400 and IMPULSE care bed with undercarriage Edition 400LR)

Use the corresponding buttons on the hand controller to adjust the bed to the low position.



(Low position operating button)

Press the button to lower the bed from the transfer position to the low position.

The undercarriages Editions 420/420LR are adjusted to the low position by means of the S8 button on the hand controller without intermediate stop.



Caution! Before pressing the button, be sure that the lowered split side rails are pushed in (☞ p. 38, Fig. 1) and there are no body parts or any objects under the bed!

Trendelenburg position

(not with IMPULSE care bed with undercarriage Edition XL and Editions 420/420LR)

Use the corresponding buttons on the hand controller to adjust the bed to the Trendelenburg position.



(Trendelenburg position button)

Care beds with HC-146/HC-147:

The Trendelenburg position is blocked for safety reasons when the bed is in resident setting.

Care staff can release the block using the key switch on the back of the hand controller (position R2 – FREE).
☞ p. 19 *Hand controller symbols*



After this, the bed can be moved into the Trendelenburg position using the S10 button on the hand controller. The incline of the sleeping surface is listed in the table on ☞ p. 14 *Technical data*.

Care beds with HB-400:

The locking function of the Trendelenburg position is activated with the aid of the magnetic chip supplied (☞ p. 22, *Locked Function Hand controller HB-400*). The locked function is released in the same way.



Caution! The Trendelenburg function cannot be implemented if a power supply failure occurs in combination with an empty battery or if the lifting motors fail. In this case, the resident must be relocated to another bed as required.



Warning! The Trendelenburg position must only be used if prescribed by a doctor. Improper use can result in lasting injury to residents.

Anti-Trendelenburg position

(not with IMPULSE care bed with undercarriage Edition XL and Editions 420/420 LR)

Use the appropriate button on the hand controller to adjust the bed to the anti-Trendelenburg position (feet lowered).



(Anti-Trendelenburg operating button)

The incline of the sleeping surface is listed in the table on
p. 14 *Technical data*.



Caution! Please consider resident safety when adjusting the bed inclination. Ensure that no body parts or any objects are located in the area of the

lifting mechanism.

Comfort sitting position

(not with IMPULSE care bed with undercarriage Editions 420/420 LR)

Use the corresponding buttons on the hand controller to adjust the comfort sitting position.



(Comfort sitting position operating buttons)

Pressing the button moves the bed quickly into a comfortable sitting position by simultaneously adjusting the sleeping surface and the lifting mechanism.

This function must only be used with mobile residents and residents without any physical problems.



Caution! Please consider resident safety when adjusting the bed to the comfort sitting position. Ensure that no body parts or any objects are located in the area of the lifting mechanism.

Hand controller locking function

The electrical unit combines state-of-the-art technology and single-fault safety.

The locking function is a further safety precaution.

Operating locking function hand controller HC-146/HC-147:

The locking function is located on the back of the hand controller and can be operated by staff with a key switch. **By turning the key switch to the various switch positions, it is possible to restrict the hand controller functions. Symbols ➤ p. 19-20.**

Operation locking function hand controller HB-400:

To activate the locking function, place the magnetic chip supplied in the area marked red. At the same time, press the respective blocking button (for secondary symbols, see table ➤ p. 22-20). A red lighted LED in the middle of the row of buttons confirms that the function has been locked.

The locked functions are released in the same way. By placing the magnetic chip and simultaneously pressing the button with the lock symbol, the function is released again and the respective LED goes out ➤ p. 22-20.

Braking and transport

IMPULSE care bed with undercarriage Edition 100

The sleeping surface is in the lowest position. Enable the castors by pressing the S8 button on the hand controller and pressing the castor button at the same time.



(Enable castors operating buttons)



The IMPULSE Edition 100 model is now in transport mode. Press the S7 button on the hand controller (height adjustment up) to retract the castors and position the bed on its feet.



(Retract castors / brake operating buttons)



On the IMPULSE care bed with undercarriage Edition 100 (only mobile in low position), the S8 button on the hand controller has a double function.



Caution! Transport mode must only be used whilst the bed is being moved. After moving the bed, disable the castors by moving the sleeping surface back to a raised position.

IMPULSE care bed with undercarriage Edition 400/420

The IMPULSE care beds with undercarriage Edition 400 and Edition 420 (mobile at any positioning height) have one central castor brake per axle, which is operated mechanically using a foot pedal accessible from both sides.



The IMPULSE care beds with undercarriage Edition 400 and Edition 420 have two different adjustment options:

1. Castors braked (foot pedal down)
2. Castors enabled for 360° movement (foot pedal up)

IMPULSE care bed with undercarriage Edition 400LR and Edition XL

The IMPULSE care beds with undercarriage Edition 400LR and Edition XL have one central castor brake, which is operated mechanically with a central foot pedal. The foot pedal is located in the central foot area of the chassis.



The IMPULSE care beds with undercarriage Edition 400LR and Edition XL have three different adjustment options:

1. Castors centrally braked (foot pedal down)
2. 4 castors enabled for 360° movement (foot pedal in the middle)
3. 1 castor set to a fixed direction (foot pedal up)

IMPULSE care bed with undercarriage Edition 420LR

The IMPULSE Edition 420LR model with optional castor sizes of 100 mm or 125 mm (mobile in any positioning height) has one central castor brake, which is operated mechanically with a central foot pedal. The foot pedal is located on both sides near the foot end axle.



The IMPULSE care beds with undercarriage Edition 420LR model has three different adjustment options:

1. Central braking of the castors (foot pedal at the bottom)
2. 4 castors enabled for 360° movement (horizontal foot pedal)
3. 1 castor set to a fixed direction (external foot pedal up)



Caution! The IMPULSE care beds with undercarriage Edition 400/400LR and Edition 420/420LR, as well as Edition XL can be moved with the sleeping surface at any height. It should only be moved with the sleeping surface raised in exceptional cases and under the supervision of care staff. Once the bed has been moved to its final position, check that the castor brakes have been applied to ensure that the bed is in braked mode. The resident's safety is the top priority!

Adjust DS side rails

Full-length side rails

In the down position, the side rails are located on top of one another next to the sleeping surface frame.

1. Raising function:

Lift the upper side rail by the side rail groove (B) until you hear the safety button (A) lock into place on both sides.




Caution! Check that the side rail is locked in place by rattling it several times.

2. Lowering function:

Slightly lift the upper side rail by the side rail groove (B), while at the same time pushing in the safety button (A). Slowly lower the side rail to the lowest position. Repeat the process at the other end of the bed.



Caution! When lowering and raising the side rails, please be extremely careful not to trap fingers, hands or any other parts of the body between the side rails and the sleeping surface frame.



Caution! For disorientated or undernourished residents, we strongly recommend using side rail padding to prevent limbs from becoming trapped between the side rail gaps, which could lead to injury.

DS side rail dimensions

Adjust GS side rails

Split side rails

1. In standby position, the side rails are located next to the sleeping surface. In this position, they prevent the mattress from slipping. (Fig. 1)
2. The side rails are raised by swivelling them upwards. In the middle position, they prevent the resident from falling out of bed. When raised at the head end, they also provide standing and mobility assistance for the resident. (Fig. 2)
3. The upright side rail is released by means of two spring catches in the side bars below the height adjustment handrail and positioned at its maximum height. Perform telescopic adjustments using both hands simultaneously when moving the rail upwards or downwards to prevent jamming. (Fig. 3) **Do not force the movement!**
4. To lower the telescoping side rail height extension, proceed the same way as to raise it.
5. To swivel the side rails back to the standby position, press the indicated latches on the lower side rail bar inwards simultaneously and initiate the tilting motion. (Fig. 4)



Caution! When raising the side rails and side rail height extensions, ensure and check that the latches engage securely. Always use both hands to move this element!



Caution! Activate the low position on the IMPULSE care bed with undercarriage Edition 400/400LR and Edition 420/420LR only with retracted side rails! (Fig. 1)



Fig. 1



Fig. 2



Fig. 3



Fig. 4

VGS side rail dimensions

Adjust VGS side rails

Vertically retractable split side rails

In standby position, the side rails are located next to the sleeping surface. In this position, they prevent the mattress from slipping.

1. To raise the side rails, pull them upwards with both hands until they engage with an audible click. In the middle position, they prevent the resident from falling out of bed. When raised at the head end, they also provide standing and mobility assistance for the resident. (Fig. A, 1.)
2. The side rail extension can be released using the two spring catches below the handrail for telescopic movement and adjusted to its maximum height. Pull the rail upwards using both hands simultaneously (Fig. B, 1./2). Take care not to tilt and jam the element. Follow these actions in reverse to lower the rail.
Do not force the movement!
3. To lower the side rails, push the two release slides inwards with both hands (Fig. A, 2) and lower the side rails carefully into the standby position.




Fig. A




Fig. B



Caution! When raising the side rails and side rail height extensions, ensure and check that the latches engage securely. Always use both hands to move this element!

Adjust EVGS side rails

Vertically retractable split side rails with one-handed operation

In standby position, the side rails are located next to the sleeping surface. In this position, they prevent the mattress from slipping.

1. To raise the side rails, pull them upwards with both hands until they engage with an audible click. In the middle position, they prevent the resident from falling out of bed. When raised at the head end, they also provide standing and mobility assistance for the resident. (Fig. A, 1.)
2. The side rail extension can be released using the two spring catches below the handrail for telescopic movement and adjusted to its maximum height. Pull the rail upwards using both hands simultaneously (Fig. B, 1./2). Take care not to tilt and jam the element. Follow these actions in reverse to lower the rail. **Do not force the movement!**
3. To lower the side rails, pull the release handle forward with one hand (Fig. A, 2.) and carefully lower the side rails into the standby position.



Caution! When raising the side rails and side rail height extensions, ensure and check that the latches engage securely. Always use both hands to move this element!



Fig. A








Fig. B

VGS side rail dimensions

When using the ABS/metal lattice sleeping surface



When using the comfort sleeping surface



Integrated bed extension

Beds equipped at the factory with an integrated bed extension (optional) can be extended by up to 20 cm without tools. For this to be possible, the side rail bars of full-length side rails must be replaced.

This function is activated in three stages and without tools using the two locking bolts at the bottom foot end of the sleeping surface:

1. On the IMPULSE care beds with undercarriage Editions 400/400 LR, pull both locking bolts upwards and then rotate them about a quarter turn to the right until they reach the locking point. On the IMPULSE care beds with undercarriage Editions 420/420 LR, the locking bolts point laterally inwards. The mode of operation is analogous to the chassis editions described above. The bed extension is now unlocked.
2. Reach below the foot board and carefully pull out the bed extension by approx. 10 cm or 20 cm.
3. Finally, turn the locking bolts back to their original positions. Afterwards, push the bed extension back carefully until the mechanism engages.



Caution! A mattress insert (accessory) must be used with 20 cm extensions. To use a mattress insert, first pull the bed extension out as far as it will go. Once the insert has been positioned, continue as described in point 3.





Caution! The sleeping surface extension must only be activated by authorised specialist staff.

IMPULSE care bed with undercarriage Editions 400/400 LR and Edition XL



IMPULSE care bed with undercarriage Editions 420/420 LR



Caution! For reasons of safety, it is essential to replace existing safety elements and attachments if the sleeping surface extension is activated on beds with full-length side rails!

MAINTENANCE

The manufacturer is only liable for the safety and reliability of the product if it is serviced regularly and used in line with the safety instructions. If any significant faults are found during maintenance work which mean the safe operation of the product cannot be guaranteed, the product must be taken out of use. Maintenance work must be carried out at least once a year.



Any faults that impair the function and safety of the care bed must be resolved before the bed is used again and must be reported to the responsible person.

Only original replacement parts from Betten Malsch GmbH are permitted to be used.



Service and maintenance tasks must not be conducted when the bed is occupied. The resident or the care staff must not conduct maintenance tasks.

Procedure

1. Visual inspection

Check the welded structures for cracked weld seams and for plastic deformation and wear. The welded structures include the undercarriage and the sleeping surface with the moving interior components. Also check that all screw joints are firmly connected.

2. Level of protection and functionality check of the side rails

During the functionality check, check whether the side rails can be easily locked in place and that no impermissible wear or deformation is visible.

The spacing requirements stated in DIN EN 60601-2-52 are used to check the level of protection. This is shown by the following diagram and table.




Fig. 1 (dimensions of a single-piece side rail)




Fig. 2 (dimensions of a split side rail)

Check whether the required spacings are also complied with under load. Check dimensions A and D using a test tool in line with DIN EN 60601-2-52. The test force for dimension A is 250 N.

Dimensions		Requirement [mm]
A	The largest dimension in at least one direction between components of the side rail/handrail in all commonly used positions.	$A \leq 120$
B	Thickness of the normally used mattress without compression as stated by the manufacturer	As stated by the manufacturer
C	Distance from the upper edge of the side rail to the upper edge of the mattress without compression with the sleeping surface in a level position.	$C \geq 220$
D	Distance between head/foot board or accessories to the side rail/handrail is parallel to the sleeping surface in a level position. This also applies to extended foot sections.	$D \leq 60$ or $D \geq 318$
E	Distance between split side rails with sleeping surface in a level position.	$E \leq 60$ or $E \geq 318$

F	The largest dimension in at least one direction of every opening below the side rail	If $D \geq 250$; $F \leq 60$ If $D \leq 60$; $F \leq 120$
G	Length of the side rail(s)	$G \geq 2/3 H$
H	Inner distance between head and foot board without any extension of these parts	No requirements

3. Functionality check of brakes

Depending on the variant, the functionality of the brakes must be checked in every position. Check whether the brakes of electrical brake systems fully retract and extend.

4. Functionality check of the lifting motors

Move all lifting motors to their end position and back again. When doing so, please observe the following points:

- Any unusual noises generated
- Synchronism of the lifting motors
- Smooth operation of the lifting motors
- Correct path of the lifting motors
- Lifting motors must switch off automatically once they reach the end position

The travel path of the lifting motors may vary depending on the model variant. If in doubt, please contact our customer service.

5. Visual inspection of mains connection cables

Check the following points on the mains connection cable:

- Visually inspect and check the function of the strain relief and kink protection
- Visually inspect insulation parts
- Visually inspect connection cables (damage, crushing)
- Visually inspect mains connection plug
- Visually inspect cable hooks

6. Visual inspection of wiring

Check the following points:

- Damage to the cables
- Correct cable routing
- Plug connections and pull-out protection firmly connected

7. Visual inspection of housing

The housings must be checked for external damage and intact seals.

8. Battery replacement

The battery must be replaced in line with specifications at least every 2 years or after every case of emergency lowering.

9. Measurements in line with DIN EN 62353

Check the electrical components of the care bed as specified in DIN EN 62353. Leakage current is measured using substitute leakage test. The limit is $\leq 500 \mu\text{A}$.

10. Visual inspection grab handle for trapeze bar

Check that the plastic components and straps show no signs of damage. The grab handles must be replaced every 5 years.

Battery replacement

The 9-V block battery/batteries must be changed at 2-year intervals and after every actuation of the electrical emergency lowering function in order to guarantee the functional capability of emergency lowering. Only use brand batteries and please dispose of old batteries in an environmentally friendly manner.

The two 9-V block batteries are inserted in the motor control unit housing below the mounting panel. The pull-out slots are secured with a screw. This must be removed beforehand.



MAINTENANCE INTERVALS

Every 2 years and after every emergency lowering

- Replace the two 9-V block batteries

Annually

- Inspection and maintenance

As required

- Lubricate mechanical components
- Replace worn components if a defect occurs.

DELIVERY AND ASSEMBLY

Betten Malsch GmbH care beds are generally delivered fully assembled, or are assembled on site by company technicians or contractual partners.

Check the delivered bed against documentation for completeness and conformity.

Any defects or damage must be pointed out to the freight company immediately and noted on the delivery document.

Signing of the delivery documents by both parties is obligatory before commissioning.

If necessary, e.g. for maintenance, simple assembly procedures can also be performed by professional authorised persons.



After maintenance and servicing work has been completed, check the functionality of electrical systems.



The manufacturer offers technician training for maintenance and servicing work on their care beds. A certificate obtained as part of this training authorises the holder to carry out technical work independently on the beds.

DISPOSAL INSTRUCTIONS

The service life of the care bed is specified as 10 years if used appropriately.

Disposal instructions

- The operator must ensure that none of the components being disposed of are infectious/contaminated.
- If the bed is scrapped, the wood, plastic and metal parts used in its construction must be separated and disposed of properly.
- If you have any questions, contact your local authority, waste-disposal company or our customer service.

Electrical component disposal

- This bed is electrically adjustable and classified as a commercially used electrical device according to WEEE Directive 2012/19/EU (implemented in Germany in the electrical equipment act).
- The electrical components used are free from banned harmful substances in line with the RoHS-II Directive 2010/65/EU.
- Replaced electrical components (drives, control units, hand controllers etc.) of these beds must be treated as electronic waste in line with the WEEE Directive and disposed of properly.

Battery disposal

- Any individual removed batteries that can no longer be used must be disposed of properly as defined by Directive 2006/66/EC (implemented in Germany in the battery legislation) and must not be discarded with domestic waste.
- For information on this matter, contact your local waste-disposal company or our Service department.

In other countries outside Germany/the EU, the relevant, applicable national requirements must be observed.

ACCESSORIES (OPTIONAL)

Trapeze bar

The trapeze bar can be inserted to the left and right of the head end in the designated mounting sockets on the sleeping surface frame. Please ensure that the bolt is properly seated in the notch provided on the mount.

Its safe working load is 75 kg.



IV drip holder

The IV drip holder can be inserted in the sockets provided on the sleeping surface frame to the left or right of the head/foot end.

The IV drip holder is only intended for attaching IV drips and not for hanging up other accessories or similar objects.

The maximum load is 8 kg (2 kg per hook).



Bedside light

The bedside light is attached to the mount provided on the sleeping surface frame in the same way as the trapeze bar.



Caution! For safety reasons, the bedside light must only be used with the care bed manufacturer's original adapter and must only be fitted by authorised specialist staff.



Follow the safety instructions in the bedside light instruction manual.

Hand controller bracket

The additional, optional hand controller bracket is used to position the hand controller within reach of the resident.



Caution! The hand controller bracket is flexible and must not be used as an aid to standing up or as a handle.

Integrated bed linen holder

The integrated bed linen holder can be easily extended by pulling at the base of the foot board. This hygienic shelf directly on the bed simplifies changing the bed linen.



TROUBLESHOOTING

Malfunction	Possible cause	Possible solution
Drives cannot be operated using the hand controller	Mains cable is not plugged in	Plug in the mains cable
	Socket not live	Check socket
	Cable plug connection not firmly connected	Check plug connections on the motor and hand controller
	Hand controller or drive is faulty	Inform the operator, specialised dealer or our customer service
	Functions locked on the hand controller	Release the function on the hand controller (☞ p. 19–20)
Motorised adjustment system is not functioning properly	There is an obstruction in the adjustment area	Check moving parts and remove any obstructions
	The safe working load has been exceeded	Reduce the load

Drives cut out after a long period of operation	The adjustment time or safe working load has been exceeded and the control unit has reacted to overheating	Allow the drive system to cool down sufficiently.
Opposite function activates when operating the hand controller buttons	Motor plugs have been mixed up	Check that the cables are connected correctly or inform your operator, specialist retailer or our customer service
The side rails can no longer be properly adjusted	The mechanism is obstructed or bent	Check all moving parts and remove any obstructions or contact our customer service
Castors do not brake or cannot be rolled	Obstructions have become trapped in the castors	Remove obstructions
	The castor system is faulty	Contact our customer service team

PRODUCT SAFETY

This product bears the CE seal and therefore meets the requirements of German and European safety standards applicable to the product.

Standard	Comment
MPDG	Medical devices implementation act
VO (EU) 2017/745	Europäische Medizinprodukteverordnung (Medical Device Regulation, MDR)
DIN EN ISO 13485	Quality management for medical products
DIN EN ISO 9001	Quality management systems
DIN EN ISO 14001	Environmental management systems
DIN EN ISO 14971	Risk analysis
DIN EN 12182	Technical aids for disabled persons
DIN EN 60601-2-52	Medical electrical equipment
DIN EN 60601-1	Medical electrical equipment
DIN EN 60601-1-2	EMC – electromagnetic compatibility
DIN EN 12530/DIN EN 12531	Medical castors
DIN EN ISO 15223-1	Symbols for labelling medical products
DIN EN 1041	Symbols and information accompanying medical products
DIN 33402-1	Human body dimensions
DIN 68861-1	Furniture surfaces
BfArM recommendation	German Federal Institute for Drugs and Medical Devices

CLEANING AND DISINFECTION

Disinfection

The care bed must be disinfected regularly and at least before every change of resident. All detergents in line with DIN EN 12720 are suitable for wipe-down disinfection of the bed. The care bed must not be disinfected in inline washing systems or using water spray. Betten Malsch GmbH recommends the care products listed by the RKI for disinfection.

The detergents used for disinfection must only be used in line with the manufacturer's instructions.



Caution! Under no circumstances use abrasives, cleaning pads or stainless-steel cleaners for cleaning. Before using any disinfectants, please consider the dosage and any potential hazards that may be caused by combining them with other substances. Remove the plug from the mains socket when disinfecting the care bed and protect the drive system from moisture.



We provide separate instructions for cleaning and disinfecting our care beds.

Care of wooden parts

Malsch care beds only use wooden parts that are finished in compliance with the DIN 68861-1A standard. The aim is to produce a comfortable design, maximum functionality and a high level of practical use. To ensure you are able to enjoy this product for as long as possible, we recommend cleaning with commercial furniture cleaning products and polishes.

Even after extremely careful selection and sorting of our wooden materials, the wood is subject to a natural ageing process. Over time, environmental influences such as air humidity, heat and UV radiation can cause changes in the colour of real wood surfaces, even when they are treated. Solid wood elements are a natural product with an individual grain and characteristics. Slight colour and shading differences within the same delivery are natural and technically unavoidable. For these reasons, relative shading and colour differences and marks due to natural growth patterns in real wood do not constitute a fault and Betten Malsch GmbH can accept no liability or warranty claims for these.

SAFE DECOMMISSIONING / STORAGE

Proceed as follows to safely decommission the bed or prepare the bed for storage:

- Disconnect the bed from the power by pulling out the mains plug.
- Activate the brake system.

Storage

- Charge the (optional) battery regularly to prevent deep discharge.
- Remove any accessories such as bed lamps, trapeze bars etc.
- Cover the care beds so that the wooden parts and the frame cannot be damaged.
- Mark the storage date clearly on the bed (so that maintenance intervals can be observed).
- Lock the hand controller.



Caution! The same conditions apply to the storage location of care beds as to the working environment (temperature, humidity, heat, etc.)



The manufacturer's transport aid must be used to transport the beds.

ELECTROMAGNETIC COMPATIBILITY (EMC)

The bed is intended for operation in the electromagnetic environments listed below. The customer or user of the bed must ensure that it is used in a suitable environment.

Guidelines and manufacturers' declarations – electromagnetic emissions

Emission measurement	Compliance	Electromagnetic environment directive
RF emissions, CISPR 11	Group 1	The bed uses HF energy exclusively for its internal functions. Therefore, it produces very low HF emissions and it is unlikely that nearby electronic devices will be adversely affected.
RF emissions, CISPR 11	Class A	The bed is suitable for use in all facilities including residential areas and those that are directly connected to a public supply network that also supplies buildings used for residential purposes.
HF emissions IEC 61000-3-2	Class A	
Voltage fluctuations and flicker IEC 61000-3-3	Compliant	

Guidelines and manufacturers' declarations – electromagnetic resistance

Phenomena	EMC basic standard or test procedure	Environments in domestic health care areas	Compliance level	Electromagnetic environment directives
Static electricity discharge	IEC 61000-4-2	±8 kV contact	±8 kV contact	The flooring material must consist of wood or concrete. If synthetic materials are used, the relative humidity must be at least 30%.
		± 2 kV; ± 4 kV; ± 8 kV; ± 15 kV air discharge	± 2 kV; ± 4 kV; ± 8 kV; ± 15 kV air discharge	
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m	10 V/m	Portable and mobile radio devices must not be used closer to the care bed and its cables than the recommended distance calculated using the equation applicable to the frequency of the transmitter.
		80 MHz to 2.7 GHz	80 MHz to 2.7 GHz	
			80% AM at 1 kHz	

Magnetic fields with energetic rated frequencies	IEC 61000-4-8	30 A/m	30 A/m	Power-frequency magnetic fields must correspond to the typical value found in the business and hospital environments
		50/60 Hz	50/60 Hz	
Electrical fast transient disturbances / bursts	IEC 61000-4-4	± 2 kV	± 2 kV	The power supply quality must correspond to that of a typical business or hospital environment.
		100 kHz repetition frequency		
Surge voltage	IEC 61000-4-5	± 0.5; ± 1 kV	± 0.5; ± 1 kV	The power supply quality must correspond to that of a typical business or hospital environment.
Cable to cable				
Surge voltage	IEC 61000-4-5	± 0.5; ± 1 kV; ± 2 kV	± 0.5; ± 1 kV; ± 2	The power supply quality must correspond to that of a typical business or hospital environment.
Cable to ground				
Conducted disturbances induced by high-frequency fields	IEC 61000-4-6	3 V	3 V	
		6 V in ISM and amateur radio frequency bands	6 V in ISM and amateur radio frequency bands	
Voltage dips	IEC 61000-4-11	0% UT; ½ cycle at 0, 45, 90, 135, 180, 225, 270 and 315 degrees	0% UT; ½ cycle at 0, 45, 90, 135, 180, 225, 270 and 315 degrees	The power supply quality must correspond to that of a typical business or hospital environment. If the user of the bed requires continued functionality even in the event of power supply disruptions, it is advisable to supply the bed with power from an uninterruptible power supply or from a battery.
		0% UT; 1 cycle and 70% UT; 25/30 cycles Single-phase: at 0 degrees	0% UT; 1 cycle and 70% UT; 25/30 cycles Single-phase: at 0 degrees	
Power failures	IEC 61000-4-11	0% UT; 250/300 cycles	0% UT; 250/300 cycles	

Interference resistance of enclosure ports to high-frequency wireless communication equipment

Test frequency [MHz]	Band [MHz]	Service	Modulation	Maximum power [W]	Distance [m]	Immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 16 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	Frequency modulation ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
145						
780						
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3,4,5; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2400 to 2570	Bluetooth; WLAN 802.11 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

WARRANTY AND SERVICE

By purchasing a care bed from Betten Malsch GmbH, you have chosen a premium, high-quality product.

Malsch care beds are covered by a 24-month warranty calculated from the date of purchase.

In the event of material or manufacturing faults occurring within the warranty period, the bed will be replaced or repaired free of charge.

This excludes faults and errors caused by inappropriate handling or external influences.

Our normal terms of business and delivery apply.



If you have any questions, please contact us on the following numbers:

Customer service

Phone: +49 (0) 6626 915-100
Fax: +49 (0) 6626 915-127

info@bettenmalsch.de
bettenmalsch.com

DECLARATION OF CONFORMITY

<p> GERMAN DESIGN AWARD WINNER 2018</p> <p>DE</p> <p>EG-Konformitätserklärung</p> <p>nach der Verordnung (EU) 2017/745 des europäischen Parlaments und des Rates vom 5. April 2017 über Medizinprodukte, Anhang IV.</p> <p>Der Hersteller: Betten Malsch GmbH Rohbergstraße 9 36208 Wildesloh-Obersuhl, Deutschland Tel: +49 (0) 6626 915-100 Fax: +49 (0) 6626 915-116</p> <p>erklärt in alleiniger Verantwortung, dass die nachfolgend benannten Produkte den grundlegenden Anforderungen und Bestimmungen der Verordnung (EU) 2017/745 für Medizinprodukte entsprechen und gemäß der Anhänge II, III und VIII (Regel 1.13) der Risikoklasse I zugeordnet werden:</p> <p>Produktnamen: Pflegebett AURA Pflegebett IMPULSE Pflegebett ATYLEEN</p> <p>Basic UDI-DI: 406584BMALSCH-PKL00002V</p> <p>Die bezeichneten Produkte wurden unter Anwendung der folgenden Richtlinien und harmonisierten Normen produziert:</p> <p>Elektrische Sicherheit: IEC 60601-1 A2:2019 Mechanische Sicherheit: IEC 60601-2-52:2009+A1:2015 EMV: IEC 60601-1-2:2014 Gebrauchstauglichkeit: IEC 60601-1-6:2019+A1:2013 IEC 62364-1:2015+COR1:2016 Risikomanagement: DIN EN ISO 14971:2020-07 Richtlinie zur Beschränkung gefährlicher Stoffe RoHS Richtlinie 2011/65/EU</p> <p>Durch die Einhaltung der Bestimmungen der Verordnung (EU) 2017/745 werden die Anforderungen zur Anbringung einer CE-Kennzeichnung erfüllt. Aufgrund der Spezifikation als Medizinprodukt Klasse I werden Produkte und Verpackung spätestens ab Mai 2020 zusätzlich mit einer UDI-Kennzeichnung versehen. Eine Konformität der Produkte und Entwicklungsdokumentation sowie des QM-Systems wird durch die Zertifizierung nach DIN EN ISO 13485:2016 bestätigt.</p> <p>Bei einer mit uns nicht abgestimmten Änderung des oben genannten Produktes verliert diese Erklärung ihre Gültigkeit.</p> <p>Wildesloh, den 23.04.2021</p> <p>Betten Malsch GmbH Gesundh. Geschäftsführer: Ralf Malsch Anspruchsbef. Geschäftsführer: HGB 1510 Erfüllungsport: Wildesloh Registernummer: 07522920992 USt-Ident.-Nr.: DE 223 627 968</p> <p>Fernsehrückmeldung: Rohbergstraße 9 D-36208 Wildesloh-Obersuhl Tel.: +49 (0) 6626 915-100 Fax: +49 (0) 6626 915-116 E-Mail: info@bettenmalsch.de</p> <p>Deutsche Bank AG BIC: BFSW33HAN IBAN: 2512 0510 0064 3553 01 BIC: BFSW33HAN IBAN: 2512 0510 0064 3553 01 BIC: BFSW33HAN IBAN: 2512 0510 0064 3553 01 BIC: BFSW33HAN IBAN: 2512 0510 0064 3553 01</p>	<p></p> <p>DN</p> <p>EC Declaration of Conformity</p> <p>in accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 concerning medical devices, Annex IV.</p> <p>The manufacturer: Betten Malsch GmbH Rohbergstraße 9 36208 Wildesloh-Obersuhl, Germany Phone: +49 (0) 6626 915-100 Fax: +49 (0) 6626 915-116</p> <p>declares under its sole responsibility that the devices named below comply with the essential requirements and provisions of Regulation (EU) 2017/745 for medical devices and are assigned to risk class I in accordance with Annexes II, III and VIII (Rule 1.13).</p> <p>Product names: Care bed AURA Care bed IMPULSE Care bed ATYLEEN</p> <p>Basic UDI-DI: 406584BMALSCH-PKL00002V</p> <p>The designated products have been produced in application of the following directives and harmonised standards:</p> <p>Electrical safety: IEC 60601-1 A2:2019 Mechanical safety: IEC 60601-2-52:2009+A1:2015 EMC: IEC 60601-1-2:2014 Usability: IEC 60601-1-6:2019+A1:2013 IEC 62364-1:2015+COR1:2016 Risk Management: DIN EN ISO 14971:2020-07 Directive on the Restriction of Hazardous Substances RoHS Directive 2011/65/EU</p> <p>By complying with the provisions of Regulation (EU) 2017/745, the requirements for affixing a CE marking are fulfilled. Due to the specification as a medical device class I, the product and packaging will additionally be provided with a UDI marking from May 2020 at the latest. Conformity of the product and development documentation as well as the QM system is confirmed by certification according to DIN EN ISO 13485:2016.</p> <p>In the event of a modification of the above-mentioned product not agreed with us, this declaration loses its validity.</p> <p>Ralf Malsch Geschäftsführer / CEO / General</p> <p>TÜV CERT TÜV CERT</p>
---	---

CERTIFICATES



CERTIFICATE

Management system as per

ISO 13485:2016
EN ISO 13485:2016+AC:2016

Evidence of conformity with the above standard(s) has been furnished
and is certified in accordance with TÜV PROFI CERT procedures for



Betten Malsch GmbH
Rohbergstraße 9
36208 Wildeck - Obersuhl / Germany

scope:

Development, manufacturing and distribution of
clinic- and health care bed systems, room furnishing

Certificate registration No. **73 108 1297** Certificate valid from 2009-07-11 to **2022-09-30**

Audit report No. 4307 3534 First certification 2017-07-11





CERTIFICATE

Management system as per

DIN EN ISO 9001:2015

Evidence of conformity with the above standard(s) has been furnished
and is certified in accordance with TÜV PROFI CERT procedures for



Betten Malsch GmbH
Rohbergstraße 9
36208 Wildeck - Obersuhl
Germany

scope:

Development, manufacturing and distribution of
clinic- and health care bed systems, room furnishing

Certificate registration No. **73 108 1297** Certificate valid from 2019-10-01 to **2022-09-30**

Audit report No. 4352 8721 First certification 2006-07-31





CERTIFICATE

Management system as per

DIN EN ISO 14001:2015

Evidence of conformity with the above standard(s) has been furnished
and is certified in accordance with TÜV PROFI CERT procedures for



Betten Malsch GmbH
Rohbergstraße 9
36208 Wildeck - Obersuhl
Germany

scope:

Development, manufacturing and distribution of
clinic- and health care bed systems, room furnishing

Certificate registration No. **73 108 1297** Certificate valid from 2019-10-01 to **2022-09-30**

Audit report No. 4352 8721 First certification 2013-10-01



Item no. 91300 100002.1
DE, Stand 05/2021, Rev. 2.1
Colours may vary

Subject to technical changes

Bakare Beds Version:
Rev 1 Initial Issue Updated: 26.4.22

Bakare Beds Ltd. | Unit 1a, Bluewater Estate, Plymouth, Devon, GB
Phone: (01752) 512227 | Fax: (01752) 511117
www.bakare.co.uk | info@bakare.co.uk

Betten Malsch GmbH | Rohbergstraße 9 | 36208 Wildeck-Obersuhl | Germany
Tel.: +49 (0) 6626 915-100 | Fax: +49 (0) 6626 915-116
bettenmalsch.com | info@bettenmalsch.de

